



**FIXED**

**Die Spacer**

- 2 Coats\*                       1 Coat
- None                                 \_\_\_\_\_ coats

**Occlusal Contact**

- Tight (Touching Opp)
- Light occlusal (1/4mm out)\*
- Out of occlusion (1/2mm out)

**Interproximal Contact**

- Light
- Medium\*
- Heavy

**Gingival Embrasure**

- Nature for both anterior & posterior\*                       Closed for both anterior & posterior
- Close for anterior, open for posterior                       Open for both anterior & posterior

**Occlusal Stain**

- |                                 |                                |                                 |                                |
|---------------------------------|--------------------------------|---------------------------------|--------------------------------|
| <b>(Brown)</b>                  |                                | <b>(Orange)</b>                 |                                |
| <input type="checkbox"/> Light* | <input type="checkbox"/> None  | <input type="checkbox"/> Light* | <input type="checkbox"/> None  |
| <input type="checkbox"/> Medium | <input type="checkbox"/> Heavy | <input type="checkbox"/> Medium | <input type="checkbox"/> Heavy |

**Gingival Stain**

- Light\*                       None
- Medium                       Heavy

**Pontic Design**

- Full Ridge\*     Modify Ridge     No Contact     Point Contact     Pontic in Socket



- Scrap pontic areas on working model\*
- Only scrap pontic areas on solid model
- Do not touch pontic areas

**If Margin Unclear**

- Contact for discussion\*
- Send back
- Do the best to trim (no guarantee)
- Require new impression

**If Occlusal Space is Needed**

- Contact for discussion\*
- Adjust opposing tooth
- Adjust prep & make reduction coping
- Make metal island
- Make metal occlusal
- Adjust prep & mark die

**Path of Insertion**

- Contact for discussion\*
- Do not adjust, make as is (no guarantee)
- Adjust prep & make reduction (no guarantee)
- Adjust & mark adjacent teeth (if problem)

**Preparation to Bulky, Undercut or Bridge not Parallel**

- Adjust adjacent and mark in red\*
- Do not adjust, make as is (no guarantee)
- Adjust and mark in red
- Contact fo discussion

**Adjacent Tooth Undercut**

- Adjust adjacent and mark in red\*
- Do not adjust, make as is (no guarantee)
- Contact for discussion

**No Bite Enclosed or provided bite feedback inadequate**

- Hand mount\*
- Make ideal (no guarantee)
- Fixed cases
- Contact for discussion

**Rx Requested Porecelain Butt Margin, but No Shoulder Margin Prepared**

- Ignore the instruction make "no show metal"\*                      Contact for discussion
- Donot adjust, just process as is (no guarantee)

**Implant Abutment**

- Adjust as needed\*
- Do not adjust, just process as is (no guarantee)
- Contact for discussion

**Crown Under a Partial (optional)**

- Will provide partial for entirety of fabrication process
- Will provide single tray full-arch impression with partial impressed, and will temporarily provide partial until martix is fabricated

**SPECIAL INSTRUCTIONS**

---



---



---



---